

Flare Mixed-Mode Column: β_2 -Agonists and Amphetamines

Landon A. Wiest,¹ David S. Jensen,² Andrew J. Miles,² Andrew E. Dadson,² Matthew R. Linford¹
Department of Chemistry and Biochemistry, Brigham Young University¹ and Diamond Analytics²
1260 South 1600 West, Orem, Utah 84058, USA

Introduction

Amphetamines are a class of chemicals that can act as stimulants, decongestants, and hallucinogens, where some compounds in this category are illicit substances.¹ All are primary or secondary amines giving these compounds higher pK_a values (9.3 – 9.8). At $pH < 9.3 - 9.8$ the molecules are increasingly protonated, generally reducing retention on C18 columns, but at $pH > 9.3 - 9.8$ they are increasingly deprotonated (neutral), which can facilitate a reversed-phase retention mechanism. Because of the lack of stability of most silica-based columns at elevated pH, there are relatively few reports of the separation of these compounds in their neutral form on silica-based C18 columns. Accordingly, it is often necessary to derivatize them.² Here we show their direct analysis at elevated pH using the Flare mixed-mode column.

Many β_2 -agonists are used to treat asthma and other pulmonary diseases by relaxing smooth muscle tissue via action on the β_2 -adrenergic receptor.³ Some of them are used illegally to increase the muscle to fat ratios in livestock.⁴ These chemicals are also amines – they are basic analytes. Similar to the amphetamines, it is advantageous to operate at elevated pH when separating them by a reversed-phase mechanism.

Separations of amphetamines and β_2 -agonists were performed at pH 12 using the Flare Mixed-Mode column from Diamond Analytics. This column is the first functionalized, carbon-based phase.⁵ As this column is diamond-based, it has stability under extreme pH conditions.⁶

Experimental

Analytes: Cimaterol, tulobuterol, mabuterol and mapenterol were purchased from Sigma-Aldrich.

(St. Louis, MO). Phenylpropanolamine and methamphetamine were obtained from Restek (Bellefonte, PA).

Sample: Analyte mixtures were created in the mobile phase as ca. 1 mg/mL solutions.

Column: Diamond Analytics Flare Mixed-Mode (2.1 x 50 mm, 4 μ m, 200Å)

System: Agilent 1290 UHPLC, binary pump, DAD, ChemStation software

Injection volume: 1.0 μ L


Temperature: 35 °C

Flow rate: 0.2 mL/min

Detection: UV/Vis Diode Array Detector (254 nm)

Needle wash: 1 min with methanol

Mobile Phase: 70:30 10 mM aqueous phosphate buffer (pH 12):acetonitrile, isocratic



Results and Discussion

Figure 1 shows baseline separation of four β_2 -agonists in under three minutes. For the latter two compounds in the separation, efficiencies exceed 50,000 N/m and tailing factors are close to 1. Figure 2 shows the baseline separation of two amphetamines and three β_2 -agonists in under three minutes. Efficiencies of the latter two compounds are again in excess of 50,000 N/m and tailing factors are ca. 1.




Figure 1. Separation of four β_2 -agonists: (1) Cimaterol, (2) Tulobuterol, (3) Mabuterol, (4) Mapenterol.

Table 1. Retention of Various β_2 -Agonists

		t_r	N/m	T_f
1.	Cimaterol	0.601	14960	1.74
2.	Tulobuterol	1.164	38060	1.40
3.	Mabuterol	1.489	52580	1.15
4.	Mapenterol	2.485	73220	1.05

Table 2. Retention of Amphetamines and β_2 -Agonists

		t_r	N/m	T_f
1.	Phenylpropanolamine	0.521	15520	2.25
2.	Methamphetamine	0.760	25060	N/A
3.	Tulobuterol	1.168	37620	1.36
4.	Mabuterol	1.488	52580	1.15
5.	Mapenterol	2.485	73400	1.06




Figure 2. Separation of β_2 -agonists and amphetamines: (1) Propanolamine, (2) Methamphetamine, (3) tulobuterol, (4) Mabuterol, (5) Mapenterol

References

1. Report to accompany H.R. 5246 19 Sep 86.
2. R. Herráez-Hernández, P. Campíns-Falcó, A. Sevillano-Cabeza, *Analytical Chemistry* 1996, 68, 734-739.
3. D. P. Tashkin, L. M. Fabbri, *Respiratory Research* 2010, 11.
4. F. Sai, M. Hong, Z. Yunfeng, C. Huijing, W. Yongning, *Journal of Agricultural and Food Chemistry* 2012, 60, 1898-1905.
5. Y. Zhang, Z. Zhang, Y. Sun, Y. Wei, *Journal of Agricultural and Food Chemistry* 2007, 55, 4949-4956.
6. G. Saini, D. S. Jensen, L. A. Wiest, M. A. Vail, A. Dadson, M. L. Lee, V. Shutthanandan, M. R. Linford, *Analytical Chemistry* 2010, 82, 4448-4456.
7. L. A. Wiest, D. S. Jensen, C.-H. Hung, R. E. Olsen, R. C. Davis, M. A. Vail, A. E. Dadson, P. N. Nesterenko, M. R. Linford, *Analytical Chemistry* 2011, 83, 5488-5501.