

## FLARE C18 Mixed-Mode Column: Separation of Ginsenosides Re and Rd by LC-MS

Hung, C-H.<sup>1</sup>, Zukowski, J.<sup>2</sup>, Jensen, D. S.<sup>2</sup>, Dadson, A. E.<sup>2</sup>, Linford, M. R.<sup>1</sup>


<sup>1</sup>Department of Chemistry and Biochemistry, Brigham Young University, Provo, UT-84602, USA

<sup>2</sup>Diamond Analytics, 1260 S 1600 W, Orem, UT, 84058, USA

### Introduction


Ginseng is a commonly used herb in Asia. In 2010, approximate 80,000 tons of it were produced in Asia (China and South Korea) and North America (Canada and the United States).[1] The worldwide market for ginseng is approximately 2 billion dollars per year. Ginsenosides or panaxosides are the active compounds in ginseng. Liquid chromatography is commonly used to assay the amounts of ginsenosides in ginseng. In this application note, ginsenoside Re and ginsenoside Rd were separated by LC-MS using the FLARE mixed-mode column.

### Analytes:




1. Ginsenoside Re  
(MW: 947.15)

2. Ginsenoside Rd  
(MW: 947.15)



**Figure 1.** Separations of ginsenoside Re (1) and ginsenoside Rd (2).



**Figure 2.** Mass spectra of ginsenoside Re (A) and ginsenoside Rd (B) from the corresponding peaks (1 and 2) in the chromatogram in Figure 1.

The particles in the FLARE C18 mixed-mode (MM) core shell column are made of carbon, nanodiamond, and an amine-based polymer, and are stable at elevated temperature and pH.[2-5] The amine-based polymer in the stationary phase provides it with mixed-mode properties, allowing it to retain analytes via

## Experimental

HPLC system: Agilent 1260 Infinity coupled with an Agilent 6230 TOF LC-MS in positive ion mode

Analytes: Ginsenoside Re and Rd (ca. 0.5 mg/mL in MeOH)

*LC Conditions:*

Column dimensions: 2.1 x 50 mm (id: 16598.16-3)

Injection volume: 6.0  $\mu$ L

Flow rate: 0.4 mL/min

Temperature: 35  $^{\circ}$ C

A: 1000 mL H<sub>2</sub>O + 1 mL TFA

B: 1000 mL ACN + 0.875 mL TFA

anionic and hydrophobic interactions with a retention mechanism that varies with pH.[2] In this application note, a 2.1 x 50 mm mixed-mode FLARE C18 column was used to separate two ginsenosides using TFA as an additive.

Gradient:

0 min 100% A 0% B

15 min 0% A 100% B

15.1 min 100% A 0% B

20 min End

*MS conditions:*

Gas temperature: 300  $^{\circ}$ C

Drying gas: 5 L/min

Nebulizer: 35 psig

Capillary voltage: 3500 V

Fragmentor: 170 V

Skimmer: 65 V

OCT: 1Rf

V<sub>pp</sub>: 750V

## References

- [1] Baeg I-H, So S-H. The World Ginseng Market and the Ginseng (Korea). *Journal of Ginseng Research*. 2013;37(1):1-7.
- [2] Hung C-H, Kazarian AA, Dadson AE, Paull B, Nesterenko P, Linford MR. Guidelines for Understanding the Retention Mechanism of the Diamond Analytics Flare Mixed-Mode Column. *Diamond Analytics*.
- [3] Hung C-H, Wiest LA, Singh B, Diwan A, Valentim MJC, Christensen JM, et al. Improved Efficiency of Reversed-Phase Carbon/Nanodiamond/Polymer Core-Shell Particles for HPLC Using Carbonized Poly(divinylbenzene) Microspheres as the Core Materials. *Journal of Separation Science*. 2013;36(24):3821-9.
- [4] Saini G, Jensen DS, Wiest LA, Vail MA, Dadson A, Lee ML, et al. Core-Shell Diamond as a Support for Solid-Phase Extraction and High-Performance Liquid Chromatography. *Analytical Chemistry*. 2010;82(11):4448-56.
- [5] Wiest LA, Jensen DS, Hung C-H, Olsen RE, Davis RC, Vail MA, et al. Pellicular Particles with Spherical Carbon Cores and Porous Nanodiamond/Polymer Shells for Reversed-Phase HPLC. *Analytical Chemistry*. 2011;83(14):5488-501.