



Column Name: FLARE WP C18
Column Dimensions: 2.1 x 100 mm
Column ID: FW36-21100
HPLC System: Agilent 1260 with Agilent 6320 TOF LC-MS
Injection Volume: 6.0 μ l
Detection: TOF MS, positive mode
Flow Rate: 0.4 ml/min
Mobile Phase: A: 0.1% TFA aqueous B: 0.0875% TFA in ACN
Gradient: Time (mins) %A %B
0.00 100 0
8.00 0 100
8.01 100 0
20.00 100 0 end
Analyte: Ferritin (MW= 474 kDa)