



Column Name: FLARE WP C18
 Column Dimensions: 4.6 x 100 mm
 Column ID: 16636.1-3 wide pore
 HPLC System: Agilent 1200
 Injection Volume: 5.0 μ l
 Detection: UV at 230 nm
 Flow Rate: 1.0 ml/min
 Mobile Phase: A: 5 ml ACN, 95 ml H₂O, 0.1 ml TFA; B: 80 ml ACN, 20 ml H₂O, 0.1 ml TFA
 Gradient:

Time (mins)	%A	%B
0.00	100	0
14.00	60	40
14.01	100	0
22.00	100	0
end		

Temperature: 35 °C
 Analyte: Hemoglobin human (MW= 67 kD)