



Column Name: FLARE WP C18
 Column Dimensions: 2.1 x 50 mm
 Column ID: 16636.1-5 wide pore
 HPLC System: Agilent 1290
 Injection Volume: 1.0 μ l
 Detection: UV at 280 nm
 Flow Rate: 1.0 ml/min
 Mobile Phase: A: 5 ml ACN, 95 ml H₂O, 2.0 ml TFA (pH 0.7);
 B: 95 ml ACN, 5 ml H₂O, 2 ml TFA
 Gradient:

Time (mins)	%A	%B
0.00	100	0
10.00	0	100
10.01	100	0
12.00	100	0 end

Temperature: 60 °C

Analyte: apo-Transferrin human and BSA

